

CUNK

Because Chemistry Matters.

Dynamic and equilibrium-based investigations of CO_2 removal from CH_4 -rich gas mixtures on zeolites

<u>C. Reichenbach¹, A. Möller¹, R. Eschrich¹, D. Klank¹, A. Brandt², K. Gleichmann²</u>

¹Quantachrome GmbH, Odelzhausen/Germany

² Chemiewerk Bad Köstritz GmbH, Bad Köstritz/Germany

Introduction: CO₂ Removal from Methane

The adsorptive separation of gas mixtures along a fixed bed adsorber is one of the most effective and thus economical separation techniques. In general, the suitability of an adsorbent for the separation of a certain gas mixture cannot be deduced from texture data only, but is determined by a multitude of complex adsorption phenomena taking place within the adsorber column. Multicomponent adsorption equilibria, sorption kinetics, heat release and the flow conditions within the

adsorber are just a few parameter that define the overall separation efficiency. The pure component sorption isotherms of CO₂ and CH₄ and the breakthrough curves of CO_2 , CH_4 and CO_2/CH_4 mixtures have been measured at different temperatures on a binder containing zeolite 13X molecular sieve (Köstrolith[®] NaMSXK) in order to investigate the CO₂ removal from methane-rich gas mixtures (relevant in natural gas purification).

Single and multi-component sorption equilibria



- High pressure isotherms have been measured for CO₂ and CH₄ with an iSorb HP volumetric sorption analyzer
- The total adsorbed amount is much higher for CO₂ than for CH₄ over the whole pressure range
- The single component isotherms suggest a high separation performance of the zeolite for CO₂-CH₄ gas mixtures

dynaSim

- The simulation tool **dynaSim** was used to calculate the mole fraction of CO₂ in the adsorbed phase with respect to the mole fraction of CO_2 in the gas phase at a total pressure of 1 bar and a temperature of 293 K
- In the same way the total loading and the partial loadings have been calculated based on multi-component SIPS



Breakthrough Curves: Measurements and Simulations

Analysis conditions

Heat profiles measured within the fixed-bed during adsorption

- Pressure: 5 bar (absolute), Gas flow: 2.5 L min⁻¹ (STP), inner diameter of the adsorber: 30 mm, Temperature: 293 K
- Gas mixture 15 % CH_4 and 5 % CO_2 in He



Breakthrough Curves

- The measured breakthrough curves (BTC's) prove the good separation performance of the zeolite under authentic process conditions
- The measured BTC's are in accordance with the predictions based on the pure component isotherms

Simulation

- The linear driving force constants k_{IDF} of CO₂ and CH₄ can be calculated by fitting a model curve to the measured BTC's



40

50

60



30

20

10

Conclusions

- The investigated zeolite reveals an outstanding separation performance in CO_2 removal from methane-rich gas mixtures
- The separation efficiency under process conditions is in perfect agreement with the predictions based on pure component isotherms
- The higher kinetic coefficient found for CH₄ has no negative impact on the separation performance
- The excellent separation is exclusively caused by the strongly preferential adsorption of CO_2 compared to CH_4
- The dynaSorb BT in combination with the simulation tool dynaSim provides a powerful set-up to investigate technical separation processes under authentic conditions

www.quantachrome.de

Quantachrome GmbH & Co. KG Rudolf-Diesel-Straße 12 85235 Odelzhausen

Tel. +49 8134 9324 0 Fax. +49 8134 9324 25 info@quantachrome.de