

<b>QUANTACHROME - Übersicht Application Notes / Applikationsberichte</b>		1 = Dispersion	e = englisch
		2 = Pulver	d = deutsch
	<a href="http://www.quantachrome.de">www.quantachrome.de</a> <a href="mailto:info@quantachrome.de">info@quantachrome.de</a>	3 = poröse Stoffe (+Katalysatoren)	
	<b>Telefon: +49 (0) 8134 93240</b>	<b>Zuordnung</b>	<b>Sprache</b>
		<b>Stoffgruppe</b>	
<b>Lfd Nr.</b>	<b>Titel</b>		
318	Aspects of the 2015 IUPAC Recommendations for Physical Adsorption Characterization	2,3	e
317	The Use of Mixture and Pure Gas Adsorption Isotherms for Further Calculations and Predictions	3	e
316	Super-ambient Pressure Gas Separation Studies Using Breakthrough System Coupled With Mass Spectrometer	3	e
315	Breakthrough Curves on Carbon Molecular Sieves – Experiment and Modeling	3	e
314	Determination of Breakthrough Curves and Effective Gas Sorption Capacities	3	e
313	Intelligent Design eliminates potential Cross-contamination in Gas Sorption Instruments	2,3	e
312	Pulse Titration on the CHEMSTAR	3	e
311	Temperature programming Capabilities of the CHEMSTAR	3	e
310	Use of Ultrasound Attenuation Spectroscopy to Determine the Size Distribution of Clay Tactoids in Aqueous Suspensions	1	e
309	Porometry, Shape Factor, Contact Angle and Tortuosity	3	e
308	Physical Adsorption Characterization of Microporous Materials: Vacuum v. Flow Outgassing	3	e
307	Optimizing Gas Sorption Parameters to Minimize B.E.T. Surface Area Analysis Time	2,3	e
306	Using the Auadyne DVS to date archaeological material	2,3	e
305	Water vapor sorption by active pharmaceutical ingredients	2,3	e
304	Ions, ion pairs and inverse micelles ind non-polar media	1	e
303	How protein and fat content influence dairy gel formation	1	e
302	Pot-life determination of epoxy resin paints in disposable measurement cells	1	e
301	Accurate monitoring of gelling temperature of carrageenan dispersions	1	e
300	Rheological measurements of corrosive and/or hazardous materials	1	e
299	High precision determination of gel point by Time-Cure-Superposition	1	e
298	Determination of optimal gelling agent concentration	1	e
297	Using the Turbiscan Stability Index (TSI) to determine the Stability Diagram of a liquid dispersion	1	e
296	Using the Turbiscan Stability Index (TSI) to determine surfactant efficiency on emulsions stabilisation	1	e
295	Stability and end-use properties of personal care products	1	e

294	Stability of heavy fuel and crude oils by monitoring aggregation and sedimentation of asphaltenes	1	e
293	Overview of various stability studies of pharmaceutical dispersions	1	e
292	Electrophoretic mobility of Poly(acrylic acid)-coated alumina particles	1	e
289	Electroacoustics of Particles Dispersed in Polymer Gel	1	e
288	Acoustic spectroscopy for colloids dispersed in a polymer gel system	1	e
287	Bestimmung von spezifischen Oberflächen verschiedener Größenordnung	3	d
286	Obtaining High Resolution, Low Pressure Isotherms with a Constant Volume/VectorDose™ Method	3	e
285	Setting & mechanical properties of cement-polymer composites	1, 2	e
284	Why Microrheology is more adapted to monitor stability of weak emulsion-gels compared to mechanical rheology?	1	e
283	Real time monitoring of emulsions stability thanks to microrheology	1	e
282	Particle size and rapid stability analyses of concentrated dispersions	1, 2	e
281	On Line characterization of particle size during grinding	1, 2	e
280	Characterization of concentrated emulsions: application to formulation and quality control	1	e
279	On Line emulsification process monitoring and control	1	e
278	Green bodies fabricated by tape casting	1, 2	e
277	Stability parameter of mineral suspension	1	e
276	Stabilisation of mullite coating	1	e
275	Stability of concentrated suspension	1	e
274	Sedimentation of bentonite	1, 2	e
273	Polystyrene/silica particles by emulsion polymerization using positively charged colloidal silica	1, 2	e
272	Study of Suspension settling	1, 2	e
271	Sedimentation behaviour study by three optical methods	1, 2	e
270	Stability and flotation behaviour of silica	1, 2	e
269	Ceramic nanoparticles from nanoemulsions	1, 2	e
268	Stability and Rheology of Polyurethane foams stabilised by a polyol/silica ge	1	e
267	Macroporous ZrO2 ceramics prepared from nanoparticles	1, 2	e
266	Viscosity and thermal conductivity of copper oxide nanofluid dispersed in ethylene glyco	1	e
265	Slurry characteristics of spray-dried ceramic powders	1, 2	e
264	Polysilsesquioxane derived ceramic foams with gradient porosity	1	e
263	Physicochemical aspects of surfactants and biosurfactants for soil remediation	1, 2	e
262	Characterization of a two-aqueous phase system containing a nonionic surfactant	1	e
261	One-step preparation method of multiple emulsions	1	e
260	Properties of water-in-oil (W/O) nano-emulsions prepared by low energy emulsification	1	e
259	Mineral interactions with reference to bioflotation nad bioflocculation	1	e
258	Low energy emulsification related to the phase behavior of the system	1	e

257	Phase separation for black liquors in wood cooking	1	e
256	Online monitoring of coagulation for water treatment	1	e
255	Light backscattering as an indirect method for detecting emulsion inversion	1	e
254	Phase transition and microstructure of emulsions	1	e
253	Studies of formation of W/O nanoemulsion	1	e
252	Dextran-induced depletion flocculation in oil-in-water emulsions	1	e
251	Serum separation and structure of depletion- and bridging-flocculated emulsions	1	e
250	Identification of optimum formulations by Turbiscan	1	e
249	Probing floc structure with different methods	1	e
248	Investigation of microstructural changes of concentrated latex sphere suspensions	1	e
247	On line study of emulsification process	1	e
246	Kinetics of droplets aggregation	1	e
245	Dependence of creaming on particle size and concentration	1	e
244	Creaming of flocculated emulsions	1	e
243	Film formation analysis by diffusive wave spectroscopy	1	e
242	Optical film formation analysis	1	e
241	Film formation analysis by optical methods	1	e
240	Film formation of coatings studied by DWS	1	e
239	New technique for the analysis of microstructure dynamics, drying processes and coating formation	1	e
238	New technology for advanced drying analysis of coatings	1	e
237	Stability of high internal phase emulsions	1	e
236	Preparation of core cross-linked micelles	1	e
235	Characterization of PVA suspension	1	e
234	Migration behaviour of Polystyrene microspheres	1, 2	e
233	Colloidal stability of Polyurethane dispersion	1	e
232	Chemical Modification of Polydienes in Latex Medium	1	e
231	Flocculation with polymer flocculants	1	e
230	Stability analysis of Non-aqueous Emulsions for elastometric films	1	e
229	A new approach for the characterization of insoluble amphiphilic copolymers	1	e
228	Stabilisation of emulsions by PVA	1	e
227	Electrical field demulsification of w/o emulsions	1	e
226	Stability of binary and ternary particle suspensions:	1	e
225	Demulsification of crude oil emulsions: a review	1	e
224	Asphaltene adsorption on minerals: Suspension stability	1	e
223	Interfacial properties of Metalworking fluids	1	e
222	Stabilisation of asphaltenes	1	e

221	Effect of water quality on the properties of drilling fluids	1	e
220	Stability of diesel-bioethanol fuels	1	e
219	The effect of emulsifier concentration on the lubricating properties of oil-in-water emulsions	1	e
218	Mechanism of crude oil demulsification	1	e
217	Lubricant emulsion stability	1	e
216	Novel method to study oil stability	1	e
215	Formulation of a bitumen emulsion	1	e
214	Role of polyelectrolyte dispersant in the settling behavior of suspensions	1	e
213	Development of new green demulsifiers for oil production	1	e
212	Stability of carbon black dispersion	1, 2	e
211	Phosphor inks for ink-jet process in PDP	1	e
210	Encapsulation of multi-walled carbon nanotubes	1	e
209	Stability of pristine carbon nanotubes	1	e
208	The preparation of hydrophobic silver nanoparticles	1, 2	e
207	Stability of nanofiber aqueous colloids	1	e
206	Stability of multi-walled carbon nanotubes	1	e
205	Transparent film of carbon nanotubes	1	e
204	Development of Nano pigment for screen display	1	e
203	Settling Characteristics of Chemical Mechanical Polishing Slurries	1, 2	e
202	Stability of nanosphere for ocular treatment	1	e
201	pH sensitive MWCNT for drug delivery	1	e
200	Stability of emulsions for parenteral feeding	1	e
199	Novel self-assembling nanogels: Stability and lyophilisation studies	1	e
198	Polysaccharide nanoparticles	1	e
197	Parenteral emulsion stability	1	e
196	Stability of Insuline W/O/W emulsion	1	e
195	Characterisation for formulation design	1	e
194	Stability of inhaler suspensions	1	e
193	Surfactant particle interaction and stability	1	e
192	Sedimentation characterization of ophtalmic suspensions	1	e
191	Following Kinetic formation of drug by backscattered light	1	e
190	Formulation of submicronic emulsions by using high pressure homogenizer	1	e
189	Study of emulsion stabilization using polymers	1	e
188	Effect of dilution on drug dispersion	1	e

187	Online particle size measurement of multiple emulsions	1	e
186	On line monitoring of emulsions for scale-up	1	e
185	Stability of pickering emulsions	1	e
184	Stability of cosmetic emulsions	1	e
183	Effect of oil content and processing conditions on the thermal behaviour and physicochemical stability of oil-in-water emulsions	1	e
182	Modified egg yolk emulsifying properties	1	e
181	Stability of acidified milk drinks	1	e
180	Effect of Xanthan and Weighting agents on stability of emulsions	1	e
179	Stability of emulsions containing sodium caseinate and dextran sulfate	1	e
178	Stability of food foams	1	e
177	Filtration of food salting solutions	1	e
176	Microencapsulated sunflower seed oil	1	e
175	Physico chemical characteristics of Recombined dairy cream	1	e
174	Phase separation of starch/xanthan mixtures	1	e
173	On line monitoring of foam	1	e
172	Electrostatics for Emulsion flocculation	1	e
171	Foam stability of milk containing dispersions	1	e
170	Sodium caseinate emulsion and foam	1	e
169	Impact of Weighting agents and sucrose on separation of beverage emulsions	1	e
168	Particle size and stability of milks	1	e
167	Phase separation on aqueous starch systems	1	e
166	Measurement of layer formation in dairy emulsion	1	e
165	Lactoglobulin effect on depletion floc	1	e
164	Phase stability of concentrated dairy products	1	e
163	On Line bubble size measurement using a multiple light scattering sensor	1	e
162	Mesostructure of fibrillar bovine serum albumin gels	1	e
161	Untersuchung des Hydratationsverhaltens von Zementen mittels akustischer Spektrometrie und Elektroakustik	1, 2	d
160	Zeta potential measurement at ceramic slurries with original concentration by means of the DT-310 / DT-1202	1, 2	e
159	Zetapotenzialmessung an keramischen, originalkonzentrierten Dispersionen mit dem DT-310 / DT-1202	1, 2	d
158	Electric conductivity of organic liquids	1	e
157	Leitfähigkeiten nichtwässriger Flüssigkeiten	1	d
156	Combined analysis of the size and shape of particles in the pharmaceutical industry	1, 2	e
155	Choosing the most suitable shape parameter in morphological analysis	1, 2	e

154	Choosing the most suitable theory in laser diffraction Mie or Fraunhofer?	1, 2	e
153	On-line measurement of the size of silicon carbide (SiC) nanoparticles	1, 2	e
152	Particle size analysis for drilling in petroleum industry	1, 2	e
151	Particle size analysis of cocoa powders	1	e
150	ISO Compliance of the dry jet dispersion technology	1, 2	e
149	Particle size analysis of emulsion in the dairy industry	1	e
148	The importance of quality methods in the precise sizing of microsphere's diameter	1, 2	e
147	Measurement of pigments with dry jet dispersion technology	2	e
146	Size and shape characterisation of plates particles for paint industries	1, 2	e
145	Size and shape analysis of metallic particles used in battery industries	1, 2	e
144	Particle size analysis of the zeolite in catalysis industry	2	e
143	Measurement of certified cement with dry jet dispersion technology	2	e
142	Measure the size of cement particles using CILAS 930e Particle Size Analyzer	1, 2	e
141	Peculiarities of live cells' interaction with micro- and nanoparticles	1, 2	e
140	Particle size determination at different concentrated Al <sub>2</sub> O <sub>3</sub> and SiO <sub>2</sub> nanopowder suspensions with the ultrasonic spectrometer DT 1200	1, 2	e
139	Monitoring nano-particles in the presence of larger particles in liquids using acoustics and TEM	1, 2	e
138	Ultrasound for characterizing cement dispersions	1	e
137	Acoustics and rheology	1	e
136	Representative On-line measurement of comminution results for nanogrinding in stirred media mill	1	e
135	Evolution of water-in-oil emulsion controlled by droplet-bulk ion exchange. Acoustic, electroacoustic, conductivity and image analysis	1	e
134	Ionic properties of so-called "non-ionic" surfactants in non-polar liquids	1	e
133	Ultrasound for characterizing liquid based food products.1. Acoustic Spectroscopy	1	e
132	Characterization of CMP slurries.Part 2. A new composite method comprised of acoustic and electroacoustic spectroscopy and sedimentation monitored with ultrasound	1	e
131	Announcement of book "Ultrasound for characterizing colloids"	1	e
130	A new way to characterize stability and performance of cosmetic emulsions and suspensions	1	e
129	Influence of chemical composition on the acoustic properties of homogenous liquids	1	e
128	Particle size distribution and micro-rheological properties of the structured concentrated dispersions	1	e
127	Surfactant titration of Kaolin slurries using $\zeta$ -Potential probe	1	e
126	Characterization of mixed dispersions by means of acoustic spectroscopy	1	e
125	Electroacoustics phenomena In concentrated dispersions	1	e
124	Acoustics and electroacoustics for emulsions	1	e
123	Acoustics and electroacoustics for ceramics	1	e

122	Electroacoustics for concentrated dispersions	1	e
121	Characterization of chemical polishing materials by means of acoustic spectroscopy	1	e
120	Sensitivity of backscattering to size variation	1	e
119	Effect of the stirring power on the emulsification efficiency	1	e
118	Characterisation of phase inversion with optical method	1	e
117	Effect of the homogenisation speed on the emulsification efficiency	1	e
116	Charge effect of flocculating agents on the flocculation efficiency	1	e
115	Effect of flocculating agent concentration on the flocculation efficiency	1	e
114	Effect of the positive ion nature on the carbonate dissolution efficiency	1	e
113	Soft matters end use properties analysis thanks to microrheology	1	e
112	How to set up a stability method	1	e
111	Kinetics characterization of two gelation processes	1	e
110	Formulation of a concentrated emulsion for natural agrochemical	1	e
109	Effect of the water hardness on dispersions	1	e
108	Stability analysis of drilling fluids thanks to microrheology	1	e
107	Stability of petroleum oil	1	e
106	Quality control of self-emulsifying oils and lubricant emulsions	1	e
105	Dispersability of multi-walled carbon nanotubes	1	e
104	Stability of platinum nanoparticles in fuel cells	1	e
103	CMP slurry stability	1	e
102	Stability of suspensions for electronic applications	1	e
101	Stability of foams	1	e
100	Foamability of a surfactant	1	e
99	Detergency properties of surfactants	1	e
98	Control of the quality and performance of surfactants	1	e
97	Effect of electrolyte introduction on injectable emulsion stability	1	e
96	Effect of an antibiotic on the stability of injectable emulsions	1	e
95	Stability of pharmaceutical products	1	e
94	Development of pharmaceutical products	1	e
93	Loading and delivery of a pharmaceutical drug in an hydrogel	1	e
92	Rheological analysis of toothpastes	1	e
91	Effect of a process on the viscoelastic properties of a body cream	1	e
90	Behaviour of a 'no-clump' mascara	1	e
89	Thixotropy of gloss lipstick	1	e
88	Behaviour of a 'no transfer' face foundation	1	e
87	Development and control of cosmetics products	1	e

86	Study of the gelling of an emulsion	1	e
85	Study of the release of an encapsulated active in a double emulsion	1	e
84	Pre-formulation of a double emulsion	1	e
83	Study of the slow breaking of a foam	1	e
82	Study of the quick breaking of a foam	1	e
81	Formulation of direct emulsions	1	e
80	Identification of the instability phenomenon of an emulsion	1	e
79	Characterization of an emulsion - quick and simple fingerprint of emulsions	1	e
78	The effect of ingredient quality on final formulation stability	1	e
77	Effect of raw materials on a formulation	1	e
76	Optimizing the cost of a low-fat emulsion	1	e
75	Optimizing the texture of a low-fat emulsion	1	e
74	Study of yogurt textural properties	1	e
73	Viscoelastic characterization of cheese thanks to microrheology	1	e
72	Development and control of beverage emulsions	1	e
71	Study of the making of sponge cake	1	e
70	Stability of whipped egg	1	e
69	Formulation of a chocolate milk	1	e
68	Study of yoghurt formation	1	e
67	Effect of the temperature on the clarification of the wort in the beer making process	1	e
66	Effect of fat content on the creaming of milk	1	e
65	Stability of various beverage emulsions	1	e
64	Drying mechanism of water-based coatings	1	e
63	Influence of the change of the pigment on the film formation	1	e
62	Finding the right open time	1	e
61	Switching to an environmentally friendly drier	1	e
60	Choosing the right coalescent	1	e
59	Curing of a UV-ink	1	e
58	Optically measuring drying of black or blue inks	1	e
57	Stability and use properties of pigment based inks	1	e
56	Study of the stability of a ceramic ink	1	e
55	Effect of the introduction of polymer on pigment dispersion stability	1	e
54	Stability of pigment inkjet inks	1	e
53	Automated software assistant for the proper calculation of BET area of microporous materials	3	e
52	Application of Quenched Solid Density Functional Theory (QSDFT) for pore size analysis of cylindrical and spherical pore carbons	3	e

51	Adsorptives for physisorption experiments: selection and their physical properties	3	e
50	Low surface area analysis by Krypton adsorption at 77.4K	2, 3	e
49	Physical adsorption characterization of mesoporous materials: Comments on hysteresis scanning, Part 1	3	e
48	Measuring pore sizes in walls of hollow fiber membranes by porometry	3	e
47	Accurate pore sizing of track-etched membranes by liquid expulsion porometry	3	e
46	Differentiating air filters by capillary flow porometry	3	e
45	Importance of zeta potential in wastewater treatment	1	e
44	Kinetics of water vapor sorption by apple	2	e
43	Water vapor sorption by milk powders	2	e
42	Application of Customized Optistat® DN Cryostat for advanced adsorption and structural characterization studies	3	e
41	Using mercury intrusion porosimetry to characterize specialty papers	3	e
40	Thermal conductivity detectors (tcd's) for flow-based measurements	3	e
39	Application of QSDFT (Quenched Solid Density Functional Theory) - A novel density functional theory for an accurate pore size analysis of disordered porous carbons	3	e
38	Aspects of a novel method for the pore size analysis of thin films based on krypton adsorption at liquid argon	3	e
37	Pycnometers aid osteoporosis and bioscaffold research	2	e
36	Hydrogen adsorption: experiment and application	3	e
35	Comments on the proper selection of adsorption or desorption branches for mesopore size analysis	3	e
34	Micropore size analysis of porous carbons using CO <sub>2</sub> adsorption at 273.15 K (0°C)	3	e
33	Some aspects of Quantachrome's NOVA (NO Void Analysis) technology	2, 3	e
32	Pore network modeling from mercury intrusion/extrusion porosimetry	3	e
31	Practical methods to prevent sample elutriation in vacuum volumetric gas sorption analyzers	2, 3	e
30	Pore size analysis by gas adsorption Part I: Aspects of the application of Density Functional Theory (DFT) and Monte Carlo simulation (MC) for micro/mesopore size analysis	3	e
29	Tapped density in the pharmaceutical industry	2	e
28	Elevated saturation vapor pressure of liquid nitrogen	2, 3	e
27	Attaching a mass spectrometer to the ChemBET 3000	3	e
26	State-of-the-art zeolite characterization: argon adsorption at 87.3 K and non local density functional theory (NLDFT)	3	e
25	Percent of solids in a slurry - determination by density measurements	1	e
24	Catalyst evaluation by temperature programmed reaction (TPR)	3	e
23	Geometric (Envelope) density without mercury	3	e
22	Advantages of using thermistor and motorized dewar for LN <sub>2</sub> coolant control	2, 3	e
21	Alternate method to determine density of cement	2	e
20	Achieving high accuracy and precision in helium pycnometry	3	e

19	Submicropore analysis	3	e
18	Supercritical gas adsorption using gas pycnometry	3	e
17	Micropore area by t-plot method	3	e
16	A method for the estimation of micropore volume and micropore surface area	3	e
15	Effect of coating of catalyst supports in mercury porosimetry	3	e
14	Sample conditioning	2, 3	e
13	Hysteresis, entrapment and wetting angle in mercury porosimetry	3	e
12	Contact angle in mercury porosimetry	3	e
11	Isotherms from mercury porosimetry	3	e
10	Pore spectra from mercury porosimetry	3	e
9	Hysteresis in mercury porosimetry	3	e
8	Desorption isotherms by continuous flow	2, 3	e
7	BET constant and site occupancy	3	e
6	Chemisorption-surface titration	3	e
5	Representative sampling for surface area measurements	2, 3	e
4	Particulate structure	2	e
3	Low surface area measurements using Krypton	2, 3	e
2	Low surface area measurements	2, 3	e
1	Small volume calibration	2, 3	e